





Heating Resistor  
Forms (Front View)



Filter:



Insulation:



Top view

(To Scale on 11x8.5" Paper)

To ensure the *Reader* is operating properly, control sets will be tested as part of a **Run Cycle**. The Software will prompt you what tubes should be loaded when and where.

Do NOT open +/- control tubes.  
Patient samples should NOT be opened once saliva sample is collected. Open tubes could contaminate future runs.

**QC Run:** Pull (●) and (+) from freezer and allow to thaw. Load (●) tube into slot 1 and (+) tube into slot 3. Press [Run] **immediately** after tubes inserted.




**Get Sample:** Take sterile loop, touch to your tongue (or right under tongue, then place it in a sample tube (unlabeled). Dispose of the loop and close the tube.

Place hinge forward so it is easier to remove tube without opening it. Never remove the middle tube (it is the temperature sensor)

**Run Cycle.** Insert the patient sample into the slot that the software directs you. Immediately insert the correct control (either ● or +) into the other slot and press Run.

**End of Day:** The software will direct you to do a final QC run. Just insert the tubes as instructed, and you can walk away, the reader will shutdown after the run.



## Quick Reference Instructions