

Supplementary Table 1. Assessment of pooled specimens with 1-12 viral copies per microliter of specimen, in pool sizes ranging from 1 to 24. To generate each pool, 200 μ L of the known positive stock solution of 1-12 copies/ μ L of specimen was combined with known negative saliva matrix. For each RNA extraction, 200 μ L of each sample was used.

Viral RNA Copy/ μ L	N	Pool size	Mean Ct values (SD; N=3)				% Triggering Reflex
			MS2	N gene	Orflab	S gene	
12 copies	3	1	28.83 (0.57)	25.4 (0.3)	26.81 (0.38)	26.01 (0.5)	N/A
	3	8	29.63 (0.24)	28.42 (0.32)	30.24 (0.33)	29.03 (0.24)	100%
	3	16	28.72 (0.34)	28.52 (0.19)	30.22 (0.1)	29.32 (0.13)	100%
	9	24	28.68 (1.28)	29.08 (0.95)	31.54 (0.66)	29.96 (1.02)	100%
	3	36	28.41 (0.23)	30.51 (0.56)	32.5 (0.69)	32 (0.67)	100%
	3	48	29.91 (0.33)	30.91 (0.42)	33.3 (0.28)	32 (0.27)	100%
6 copies	3	1	28.24 (0.45)	25.85 (0.32)	27.74 (0.18)	27.74 (0.77)	N/A
	3	8	28.11 (0.21)	29.31 (0.68)	31.1 (0.65)	30.11 (0.8)	100%
	3	16	28.05 (0.44)	29.45 (0.38)	31.35 (0.24)	29.95 (0.64)	100%
	9	24	29.1 (0.6)	30.39 (0.68)	32.87 (0.97)	31.78 (0.68)	100%
	3	36	28.59 (0.46)	31.96 (0.28)	34.23 (0.34)	32.75 (0.72)	100%
	3	48	28.84 (0.34)	31.79 (0.61)	33.56 (0.97)	31.59 (0.29)	100%
3 copies	3	1	27.54 (0.66)	27.15 (0.94)	28.92 (0.73)	27.74 (0.88)	N/A
	3	8	27.35 (0.44)	29.21 (0.87)	32.74 (0.23)	30.88 (2.43)	100%
	3	16	27.83 (0.08)	30.5 (1.24)	32.97 (0.02)	32.08 (0.71)	100%
	9	24	28.17 (0.41)	30.68 (0.47)	32.99 (0.99)	32.08 (0.8)	100%
	3	36	28.88 (0.34)	31.82 (0.22)	34.18 (NA)	31.72 (NA)	100%
	3	48	28.02 (0.3)	31.1 (0.18)	34.38 (0.21)	33.58 (0.27)	100%
1 copy	3	1	28.81 (0.63)	28.2 (0.1)	30.2 (0.31)	29.5 (0.34)	N/A
	3	8	28.87 (0.69)	30.15 (1.17)	33.3 (0.57)	31.24 (0.73)	100%
	3	16	28.36 (0.7)	29.52 (0.76)	33.67 (0.9)	31.83 (0.62)	100%
	9	24	27.55 (0.61)	31.11 (0.29)	33.19 (0.51)	32.6 (0.38)	100%
	3	36	28.73 (0.62)	31.72 (0.81)	34.11 (0.09)	32.12 (1.46)	100%
	3	48	28.61 (0.47)	31.61 (0.61)	34.22 (0.18)	33.62 (2.05)	100%