

Supplementary Information:

Serological Response in Lung Transplant Recipients after Two Doses of SARS-CoV-2 mRNA Vaccines

Madhusudhanan Narasimhan, PhD, MSHA^{1#}, Lenin Mahimainathan, PhD, MB(ASCP)^{1#}, Andrew E. Clark, PhD, DABMM¹, Amena Usmani, MD¹, Jing Cao, PhD, DABCC¹, Ellen Araj, MD¹, Fernando Torres, MD², Ravi Sarode, MD^{1,2}, Vaidehi Kaza, MD, MPH², Chantale Lacelle, PhD, D(ABHI)^{1*}, Alagarraju Muthukumar PhD, DABCC^{1*}

¹ Department of Pathology, University of Texas Southwestern Medical Center, Dallas, Texas

² Department of Internal Medicine, University of Texas Southwestern Medical Center, Dallas, Texas

- Equal Contribution

***Corresponding author:**

Alagarraju Muthukumar, PhD, DABCC or Chantale Lacelle, PhD, DABHI

UT Southwestern Medical Center, Dallas

Email: alagarraju.muthukumar@utsouthwestern.edu (AM) (or)

chantale.lacelle@utsouthwestern.edu (CL)

Ph: 214-645-5103 (AM)/214-648-0904 (CL)

Fax: 214-645-5106




Figure S1

Table S1: IgM_{SP}, IgG_{NC}, and IgG_{SP} Responses following Double Dose COVID-19 Vaccine Regimen in Lung Transplant Patients, Stratified by Cylex ImmuKnow assay levels

Information	n (%)	Cylex ImmuKnow assay levels		
		Low n (%)	Moderate n (%)	Strong n (%)
Two dose vaccine administered lung transplanted total subjects that had Cylex results recorded	56/73 (77)	22/56 (39.3)	26/56 (46.4)	8/56 (14.3)
IgG_{SP} serology results				
BNT162b2 (Pfizer-BioNTech)	+	3/15 (20)	3/17 (18)	0/6 (–)
	–	12/15 (80)	14/17 (82)	6/6 (100)
mRNA-1273 (Moderna)	+	1/7 (14)	4/9 (44)	1/2 (50)
	–	6/7 (86)	5/9 (56)	1/2 (50)
IgM_{SP} serology results				
BNT162b2 (Pfizer-BioNTech)	+	1/15 (7)	1/17 (6)	0/6 (–)
	–	14/15 (93)	16/17 (94)	6/6 (100)
mRNA-1273 (Moderna)	+	0/7 (–)	0/9 (–)	0/2 (–)
	–	7/7 (100)	9/9 (100)	2/2 (100)
IgG_{NC} serology results				
BNT162b2 (Pfizer-BioNTech)	+	0/15 (–)	1/17 (6)	0/6 (–)
	–	15/15 (100)	16/17 (94)	6/6 (100)
mRNA-1273 (Moderna)	+	0/7 (–)	0/9 (–)	0/2 (–)
	–	7/7 (100)	9/9 (100)	2/2 (100)

+ and – , Serology assay's positive and negative results based on the manufacturer recommended corresponding cut-off value.