

***Plasmodium falciparum* hrp2 and hrp3 gene deletion status in Africa and South America by highly sensitive and specific digital PCR**

Claudia A. Vera-Arias *et al.*

Supplementary File S1

ddPCR assays

Supermix from BioRad
All primers at 10 μ M

***hrp2* exon 2 ddPCR**

	[μ L]	Cycling conditions
SuperMix for Probes (no dUTP)	11	
hrp2_Exon 2_fwd	1.76	95° 10 min
hrp2_Exon 2_rev	1.76	94° 30 sec
hrp2_Exon 2_Probe	0.88	56° 1 min => 44 x back to step 2
tRNA_fwd	0.44	98° 10 min
tRNA_rev	0.44	4° hold
tRNA_Probe	0.22	
H2O	3.5	
DNA	2	

***hrp2* exon 1 ddPCR**

	[μ L]	Cycling conditions
SuperMix for Probes (no dUTP)	11	
hrp2-exon1_fwd	1.76	95° 10 min
hrp2-exon1_rev	1.76	94° 30 sec
hrp2-exon1_Probe	0.88	54° 1 min => 44 x back to step 2
tRNA_fwd	0.44	98° 10 min
tRNA_rev	0.44	4° hold
tRNA_Probe	0.22	
H2O		
DNA	2	

***hrp3* ddPCR**

	[μ L]	Cycling conditions
SuperMix for Probes (no dUTP)	11	
hrp3_fwd	0.66	95° 10 min
hrp3_rev	0.66	94° 30 sec
hrp3_Probe	0.44	55° 1 min => 44 x back to step 2
tRNA_fwd	0.44	4° hold
tRNA_rev	0.44	
tRNA_Probe	0.22	
H2O	3.24	
DNA	6	

Nested PCR assays

hrp2 nested PCR

Polymerase, buffer, and dNTPs from Solis Biodyne
All primers at 10 μ M

Primary PCR	[μ L]	Cycling conditions
Taq FirePol	0.15	
Buffer B	1.5	94° 10 min
MgCl ₂	1.5	94° 50 sec
dNTPs	1.3	55° 50 sec => 39x back to step 2
Primer <i>hrp2_F1</i>	0.2	70° 1 min
Primer <i>hrp2_R1</i>	0.2	70° 5 min
H ₂ O	8.25	
DNA	2	

Nested PCR	[μ L]	Cycling conditions
Taq FirePol	0.15	
Buffer B	1.5	94° 10 min
MgCl ₂	1.5	94° 50 sec
dNTPs	1.3	55° 50 sec => 39x back to step 2
Primer <i>hrp2_F2</i>	0.2	70° 1 min
Primer <i>hrp2_R1</i>	0.2	70° 5 min
H ₂ O	8.25	
DNA	2	

Primer sequences

Hrp2 F1 5'-CAAAAGGACTTAATTTAAATAAGAG-3'
Hrp2 R1 5'-AATAAATTTAATGGCGTAGGCA-3'
Hrp2 F2 5'-ATTATTACAGAACTCAAGCAC-3'

msp2 nested PCR

Polymerase, buffer, and dNTPs from Solis Biodyne
All primers at 10 μ M

Primary PCR	[μ L]	Cycling conditions
ddH ₂ O	7.85	
Buffer B	1.5	95° 1 min
MgCl ₂ 25 mM	1.5	95° 30 sec
dNTPs 2.5 mM	1.2	55° 30 sec
Primer <i>msp2_S2-fw</i>	0.4	72° 1 min => 34 x back to step 2
Primer <i>msp2_S3-rev</i>	0.4	72° 5 min
Taq FirePol	0.15	
DNA	2	

Nested PCR	[μ L]	Cycling conditions
ddH ₂ O	8.975	
Buffer B	1.5	95° 1 min
MgCl ₂ 25 mM	1.5	95° 30 sec
dNTPs 2.5 mM	1.2	55° 30 sec
Primer M5_Fc27 (10 μ M)	0.15	72° 1 min => 34 x back to step 2
Primer N5_3D7 (10 μ M)	0.225	72° 5 min

Primer S1-tail fw	0.3
Taq FirePol	0.15
Primary PCR as template	1

Primer sequences

<i>msp2_S2</i> -fw:	5' - GAAGGTAATTAAAACATTGTC-3'
<i>msp2_S3</i> -rev:	5' - GAGGGATGTTGCTGCTCCACAG -3'
<i>msp2_S1Tail</i> -fw :	5' - GCTTATAATATGAGTATAAGGAGAA -3'
<i>msp2_FC27 type_M5</i> -rev:	5' - GCATTGCCAGAACTTGAA -3'
<i>msp2_3D7 type_N5</i> -rev:	5' - VIC-CTGAAGAGGTACTGGTAGA -3'