

Supplemental Table S21. Reagents

Method	Reagent / Resource	Company	Identifier
Cell culture	BEGM Bullet kit Medium	Lonza	CC-3170
Cell culture	NHBE-Bronchial Epi Cells	LONZA	CC-2540
Cell culture	Hexa-Peptide von Adam (Tyr-Leu-Tyr....)	Thermo	PEP98UNMOD
Cell culture	HOE 140 (B2R antagonist)	SigmaAldrich	H157-250ug
Cell culture	100mg Hydrocortison	sigma aldrich	H0396-100mg
Cell culture	antibiotic antimycotic	Thermo	15240062
ELISA	human ACE2 ELISA Kit	abcam	ab23564g
ELISA	Coomassie Plus Bradford Assay Kit	thermo fisher	23236
RNA extraction	Rneasy Micro Kit (50)	Qiagen	74004
RNA extraction	QIAshredder (250)	Qiagen	79656
RNA extraction	Rnase-Free Dnase Set	Qiagen	79256
Bioanalyzer	agilent rna 6000 nano reagent	agilent	5067-1512
Microarray	Gene Expression Hybridisierungskit	Agilent	5188-5242
Microarray	Genexpression Wasch-Paket	Agilent	5188-5327
Microarray	Sure Print G3 Human GE V2 8x60K Kit	Agilent	G4851 B
Microarray	LowInput QuickAmp Labeling Kit One-Color	Agilent	5190-2305
Microarray	RNA Spike In Kit-One Color	Agilent	5188-5282
Microarray	Gasketslides	Agilent	G2534-60015
Microarray	Ethanol for molecular biology	vwr/Merck	1.085.430.250
cDNA Transcription	high capacity cDNA kit	life technologies	4368814
PCR	FS universal SYBR Green Master Mix	Sigma	4913850001
PCR	<i>Human TMPRSS2</i> fwd AGAGTGCGACTCCTCAGGTA	Metabion	
PCR	<i>Human TMPRSS2</i> rev AGGCGAACACACCGATTCTC	Metabion	
PCR	<i>Human ACE2</i> fwd ATGAAGGCCCTCTGCACAAA	Metabion	
PCR	<i>Human ACE2</i> rev CCAAGCCTCAGCATATTGAACAG	Metabion	
PCR	<i>Human ACTIN</i> fwd ATTGCCGACAGGATGCAGAA	Metabion	
PCR	<i>Human ACTIN</i> rev TCTTCATTGTGCTGGGTGCC	Metabion	
PCR	<i>Human HPRT</i> fwd TGACACTGGCAAACAATGCA	Metabion	
PCR	<i>Human HPRT</i> rev GGTCCTTTTCCACCAGCAAGCT	Metabion	
PCR	<i>Human 18S</i> fwd	Qiagen	249900
	Hs_RRN18S_1_SG QuantiTect Primer Assay		QT00199367
PCR	<i>Human 18S</i> rev	Qiagen	249900
	Hs_RRN18S_1_SG QuantiTect Primer Assay		QT00199367
PCR	<i>SARS-CoV-2 RdRP</i> fwd	Metabion	
	CGTCTGCGGTATGTGAAAG		
PCR	<i>SARS-CoV-2 RdRP</i> rev	Metabion	
	TAAGACGGGCTGCACTTACA		
PCR	<i>SARS-CoV-2 N1</i> fwd: GACCCAAAATCAGCGAAAT	Metabion	
PCR	<i>SARS-CoV-2 N1</i> rev: TCTGGTACTGCCAGTTGAATCTC	Metabion	
PCR	<i>SARS-CoV-2 N1</i> probe: 5'-FAM-ACC CCG CAT TAC GT	Metabion	