

Formula	Within 0.25 D	Within 0.50 D	Within 0.75 D	Within 1.00 D
Vs. Barrett	1.18e-1	<u>4.16e-2</u>	8.33e-2	<u>1.24e-2</u>
Vs. Haigis	<u>2.06e-06</u>	<u>9.98e-06</u>	<u>4.83e-06</u>	<u>1.84e-2</u>
Vs. Hoffer Q	<u>3.70e-11</u>	<u>3.37e-11</u>	<u>2.25e-10</u>	<u>3.96e-05</u>
Vs. Holladay 1	<u>4.83e-07</u>	<u>8.12e-07</u>	<u>9.63e-08</u>	<u>2.56e-4</u>
Vs. SRK/T	<u>1.89e-06</u>	<u>9.43e-08</u>	<u>7.98e-07</u>	<u>5.74e-06</u>

Table S2 The p-values from Cochran's Q test for the comparison of number cases within 0.25 D, 0.50 D, 0.75 D, and 1.00 D between our method and conventional methods in the testing set. All significant p-values (p-value < 0.05) were underscored. D: diopter.