

SUPPLEMENTARY MATERIAL

**The Coronavirus Calendar (CoronaCal): a Simplified SARS-CoV-2 Test System
for Sampling and Retrospective Analysis**

David S. Thaler,^{1, 2, 3, *} Manija A. Kazmi,^{3, *} Karina C. Åberg,^{3, *} Jordan M. Mattheisen,^{3,4} Thomas Huber,³ and Thomas P. Sakmar³

¹Biozentrum, University of Basel, CH-4056 Basel, Switzerland

²Program for the Human Environment, The Rockefeller University, New York, NY, USA

³Laboratory of Chemical Biology and Signal Transduction, The Rockefeller University, New York, NY, USA

⁴Tri-Institutional Program in Chemical Biology, New York, NY USA

* Authors contributed equally

Corresponding authors: E-mails: david.thaler@unibas.ch, sakmar@rockefeller.edu


Keywords:

SARS-CoV-2, COVID-19 saliva test, virus, epidemiology, public health, sample preservation

Running title:

SARS-CoV-2 sample preservation on paper diaries

Supplemental Figure 1. **The CoronaCal diary.** *Upper panel.* CoronaCal collection sheets were printed onto ordinary printer paper. Each 8-1/2" x 11" sheet can accommodate two-weeks of daily samples, and each participant was provided with enough sheets for up to 12 weeks of sampling. *Lower panel.* Covering stickers with a symptom check list were printed on Avery sticker stock. A later version of the symptoms list included "can't smell" and "can't taste." The collection sheet and sticker templates are reduced here from their actual size. Downloadable CoronaCal templates are available www.CoronaCal.com as a supplement.



Supplementary Figure 2. **RT-qPCR sensitivity plot.** As described in Materials and Methods the limit of detection and linearity of the SARS-CoV-2 virus assay was tested using the control plasmids IDT 2019 nCoV_N_Positive Control and IDT 2019 Hs_RPP30_Positive Control. Plasmids were serially diluted from 100,000 copies/rxn down to 0.1 copies/rxn in water and RT-qPCR reactions were set up in triplicate and performed as described to determine the Cq values for each dilution.