

3 **Supplemental Table 2. Reagents for ELISA Testing*****Standard Curves**

Recombinant Human Protein	Source and Catalog Number	Curve Range
C-Reactive Protein (CRP)	Abcam ab167710	15.625 to 1000 pg/ml
Calprotectin or S100A8/A9 Heterodimer	Biolegend 753404	3.9 to 250 ng/ml
High Mobility Group Box-1 (HMGB1)	Sigma-Aldrich H4652	3.9 to 250 ng/ml

4

Capture Antibody

Source and Target Species of IgG	Source and Catalog Number	Concentration Used
Rabbit polyclonal anti-human CRP	Abcam ab31156	2.5 µg/ml
Mouse monoclonal anti-human Calprotectin	LS-Bio LS-C96223	1.5 µg/ml
Rabbit polyclonal anti-human HMGB1	Upstate Biotechnology 07-584	1.5 µg/ml

5

Detection Antibody

Species	Source and Catalog Number	Concentration or Dilution Used
Mouse monoclonal IgG anti-human CRP	Abcam ab136176	1.5 µg/ml
Rabbit polyclonal anti-human Calprotectin	LS-Bio LS122793	1:1500
Mouse monoclonal anti-human HMGB1	R&D Systems MAB1690	0.75 µg/ml FC

6

Enzyme Linked Antibody

Species	Source and Catalog Number	Dilution Used
Goat Anti-Mouse Polyclonal IgG-FC	Millipore AP127P	1:2000
Goat Anti-Rabbit IgG-FC	Santa Cruz Biotechnology SC-2004	1:2000
Goat Anti-Mouse Polyclonal IgG-FC	Millipore AP127P	1:2000

7

* 96-well plates (Costar, Corning Inc, Costar, NY, USA) were incubated overnight with capture antibody diluted according to manufacturer recommendations or 1% bovine serum albumin (BSA, MilliporeSigma, Burlington, MA, USA) or 10% newborn calf serum in phosphate buffered saline (NCS-PBS, MilliporeSigma). Plates were washed 4 times with PBS prior to addition of 1st antibodies then washed 6 times prior to addition of horse radish peroxidase-conjugated 2nd antibodies. Standards were diluted 1:2 for standard curves using manufacturer's recommendations or 1% BSA or 10% NCS-PBS (MilliporeSigma). All assays were incubated with 3,3',5,5'-Tetramethylbenzidine (TMB) substrate solution (Thermo Scientific, Waltham, MA, USA). Reactions were stopped using 0.18 M H₂SO₄ after 15-30 minutes and read in an ELISA plate reader at OD₄₅₀. Standard curves were constructed using linear regression of log-transformed mean fluorescence intensities and log-transformed known protein concentrations.

15