

**Rapid Research in Diagnostics Development for TB Network
(R2D2 TB Network) Study**

Study Title: Rapid Research in Diagnostics Development for TB Network (R2D2 TB Network) Study

SOP Title: Cough sound recording

Division: SPECIMEN COLLECTION		SOP ID#: Cough sound recording
Effective date: March 2022 SOP review cycle: 6 months SOP reviewers:		Written by: Devan Jaganath, Sophie Huddart and HyfeApp team
Date Reviewed:	Approved by:	Approved Signature:

Training Implications:

- Study Coordinators
- Research Assistants
- Study Nurses

The following study personnel have read and understood the cough sound recording SOP:

Date	Name	Signature

A. Purpose

This document describes the collection of cough sounds at baseline using the Hyfe Research application on a smartphone. Hyfe Research collects and stores cough-like sounds continuously for research purposes. Procedures for quality assurance are also described.

B. Personnel needed

1. Study nurse/Research assistants

C. Supplies

1. Study-issued phone with charger
2. Carrying bag or other holder for phone that has been approved by field testing with Hyfe
3. Google Play log-in information for the study site
4. A log that keeps a record of the google log-in information, each phone, the phone password and the App Lock password
5. Spreadsheet with list of participant IDs and Hyfe PIN Numbers
6. Tablet with REDCap mobile application, including forms for Lung and Cough Sound Recording, Check-out Phone, Check-in Phone
7. Participant contact information, including a locator form
8. Tripod that can hold a phone
9. Access to internet, either via wifi or with data from the phone
10. Disposable gloves
11. Sanitizers: 70% alcohol or small batches of bleach made daily
12. Participant’s appointment card (given at enrollment)

D. Hyfe ID and Pin

1. A Hyfe ID and Pin is needed to log in to the Hyfe Research application
2. The Hyfe ID will be the same as the participant’s study ID (i.e., R2D2-country code-ID number)
3. The PIN number corresponding with the Hyfe ID will be found on the spreadsheet, for example:

Hyfe ID	Pin
R2D204001	112233
R2D204002	998877
R2D204003	565656

E. Procedures for phone set-up

**Note: this only needs to be done once for each phone*

1. Turn on the phone and use the corresponding password from the phone log


2. Delete any applications (such as games, Facebook, or Whatsapp) that you do not want the participant to use

Downloading Hyfe Research Application

1. Go to the Google Play store
2. Log in with the study site google account, specific for each study site:

Uganda	r2d2.uganda@gmail.com
Philippines	r2d2.philippines@gmail.com
India	india.r2d2@gmail.com
Vietnam	r2d2.vietnam@gmail.com
South Africa	r2d2.southafrica@gmail.com

3. Type in the search field “Hyfe Research” and click the Hyfe Research application from Hyfe
4. Press “Install” to download the application




Downloading and installing AppLock

5. In Google Play, type into the search field “AppLock”



AppLock
DoMobile Lab
In-app purchases

6. Press the “AppLock” application from DoMobile Lab and download
7. Press “Open” to open AppLock
8. AppLock will ask you to create an unlock pattern; use the pattern noted on the phone log




9. Draw the same pattern again to confirm
10. Press “Agree and Start” to the Private Policy and User Agreement
11. Once in the application, you can click the “lock” on the right of each application to block use of an application
 - a. Press the lock for the Google Play Store

Advanced



Google Play Store

Prevent installing or uninstalling application



- b. You will be asked to permit usage of AppLock, press “Permit” and then you will be taken to a screen to give permission to AppLock
 - c. You will need to reopen the application (Press “Open”) and then re-enter the AppLock pattern
 - d. Then press the lock for the Google Play Store and it should lock
 - e. Then press the lock Settings, Google Play Games, Chrome, Files, Calendar, Contacts, and the Camera
12. Exit the application

13. Now if any of these applications are opened, the AppLock screen will appear and prevent use unless the AppLock pattern is known
 - a. Locking Settings prevents deletion of the Hyfe application and doing a factory reset

F. Location and order of Baseline Cough Recording.

1. **Location:** This should be done in a well-ventilated area or negative pressure room away from other individuals. Ideally the participant should be able to remove their mask to record the baseline cough sounds.
2. **Order:** The five solicited cough sounds can be recorded at any time during enrollment, but NOT immediately (~within 2 hours) after an induced sputum as the cough sound may change. A suggested time would be before expectorated or induced sputum is being collected, as they will already be in a ventilated area.

G. Procedures with participant prior to recording

1. Don appropriate PPE, including outerwear (a lab coat, surgical gown, or apron), clean, disposable gloves, goggles, and a N-95 mask.
2. Ensure the phone is fully charged, sim card installed (if available)
3. For the baseline recording and calibration, the phone should be connected to the internet via wifi or the phone's data
4. Turn on the phone and put in the phone passcode
5. Go to settings→Apps→ Hyfe Research and press "Permissions"
6. Ensure "Microphone" is turned on and turn on "Location" if the participant consented to having their location recorded
- 7.
8. Wipe the phone and any potentially contaminated surfaces using 70% ethanol

H. Procedures for recording baseline cough sounds


1. Place study phone (turned on with Hyfe Research application running) on tripod:



2. Request the patient put his/her phone on silent mode.
3. Ensure the microphone on the phone is facing the patient



4. Let the patient stand 2-3 feet away from the phone
 - a. If the patient cannot stand, place a chair
5. Adjust the tripod height, so the phone is approximately at the same height as the patient's mouth
6. Stand at least 6 feet (2 meters) away from the patient. Angle the phone's microphone at a 45-degree angle to the patient's mouth to ensure clear recording.
7. Use a script explaining the procedure to the patient:
 - a. Voluntary coughs: *"We will ask you to cough 5 times into the microphone. Please remove your mask or shield while coughing if possible. Please try to cough as naturally as possible and try not to cover your face with your fist. You may hold onto a table or sit if that is helpful. Should you feel uncomfortable for any reason during this portion of the study, please let us know and we will work to remedy the situation."*
8. Turn on the phone and put in the phone passcode
9. Open the Hyfe Research application
10. Press "Start now"
11. Read the Disclaimer to the participant and Press "Understood"
12. Read the User Agreement and Press "Agree and Continue"
13. Sign in with the participant ID for Hyfe ID and corresponding Pin listed on the spreadsheet
 - b. Leave Optional ID blank



14. Press “Login”

The next step is to calibrate the application for the participant’s cough, which will also collect the baseline cough sounds. The phone needs to be connected to the internet for this part. This should be in a well-ventilated or open area, and continue to stay at least 6 feet (2 meters) from them while guiding them through the process.

15. Press “OK, I’m ready”

16. If asked, allow Hyfe Research to record audio

17. Press the Hyfe logo on the bottom of the screen




18. When the application says “Cough Once” the participant should cough only once directed at the microphone and wait for “Calibration” to finish.

19. Repeat the above four more times.

c. If it does not detect a cough, it will need to be completed again. Reasons for not detecting a cough include:

- i. High levels of background noise
- ii. Coughing multiple times instead of once
- iii. Coughing too early or too late

20. Once complete, press “Continue.” The application will be in continuous recording mode. It will continue to record in the background until you log out, using the same Pin as above.



21. Complete the HyfeApp section of the “Lung and cough sound recording” form in REDCap
22. Remove the phone from the tripod
23. Wipe the phone and any potentially contaminated surfaces using 70% ethanol
24. Disinfect the patient’s hands using hand sanitizer or ask the patient to wash his or her hands with soap.
25. If the participant consents to longitudinal recording, move to the next step for longitudinal recording. If they do not consent, press “logout” and enter the pin.

I. Procedures for longitudinal cough sound recording

1. Eligibility: Participants 1-300 at each clinical site will be eligible for longitudinal cough sound recording, except for Uganda where participants 1-600 will be eligible.
2. Review with the patient that the phone should remain on and near the participant at all times, but that it will only record cough-like sounds (less than a second) and not their conversations
3. The phone does not need to be connected to the internet for longitudinal recording, but will need to be connected to the internet when the phone is returned to upload the cough sounds.
4. Give the phone to the participant with the carrying device (i.e. neck bag), and guide them in turning on the phone and showing that the phone will record in the background
5. The phone should have a label indicating the phone’s ID number, enter that into the REDCap “Check-Out” form
6. Review phone use, maintenance and charging (phone should be charged daily)

7. The best place to charge the phone is in the evening on a bedside cabinet that is within 1.5 meters of the participant
8. Explain that the phone has certain applications blocked, including the ability to download new applications, use the camera or browse the internet
9. Complete the “Check out” form on REDCap
10. Provide the participant the phone, charger, and completed appointment card and remind patients to return to the clinic in two weeks to return the phone and charger
11. Write on the appointment card the Hyfe ID, Hyfe Pin and phone password in case the phone turns off and the application needs to be turned back on
12. Provide the participant the study phone number if they have any questions or problems
13. Remind participants to contact the study should a phone be damaged, lost, stolen or have lose any functionality as soon as possible.
14. Remind the participant of the importance of this data for constant recordings for two week period and that the phones are to be returned and reset so that other participants can have their data recorded.

J. Carrying and/or wearing the phone

Key points

- Phone should remain close to the participant at all times, within 1.5 meters
- The microphone should be facing up towards the participant
- There are two options that appear to work best for recording:

1. **Armband option:** If the participant chooses to wear the phone on an armband, be sure to orient the microphone up, pointing towards the mouth, like in the red circle pictured below.



2. **Necklace/bag option:** If the participant chooses to wear the phone on as a necklace, or a bag around the neck, be sure to orient the microphone up, pointing towards the mouth, like in the red circle pictured below.



K. Participant use of own phone for longitudinal collection

1. Having the participant use their own adds some challenges as there can be variability in sound recording across brands, and the Hyfe will use the participant's data to upload if their data is on or they are connected to Wifi
2. If the individual will only participate if they use their own phone, then explain to them that they can use it under the following conditions:
 - It will need to be placed in the carrying device as above
 - The cellular data should be turned off, or sim card taken out to prevent using the individual's data to upload cough sounds
 - Wifi will be turned off to prevent using data to upload cough sounds

L. Procedures for Two-Week visit

1. Ask the participant to give you the phone and charger
2. Confirm the phone is still on, functioning, and the Hyfe Research application is still running
3. Log into the phone, stop recording on the application, ensure cough data has been uploaded
 - a. If the phone has no data, either load more data or take the phone to a place with Wifi to ensure the cough sounds can be uploaded
 - b.
4. Delete any contacts, pictures or other files the participant may have stored on the phone
5. Turn-off phone and sanitize
6. Complete the "Check in" form on REDCap

M. Retention Plan

1. Reassure participants that the Hyfe Research application will not record their conversations, only short bursts of sound (less than a second) that may be a cough
2. Offer for them to keep the phone carrying bag after the study as an incentive
3. Reducing phone loss/damage

- a. The phone will be deactivated from downloading applications, using the camera or browsing the internet
- b. The participants will be reminded that the phone can be tracked if it is lost and deactivated remotely
- c. If possible, a staff member to accompany the participants to their residences to ascertain their place of residence.
4. A staff member will call the participant daily for the first three days of enrollment to discuss any issues with the phone or Hyfe Research application
5. Obtain the participant's home location using a locator form
 - a. Home visits are helpful if feasible (at least 2 times)
 - b. Interact with the participant, elders and community leaders to discuss activities and importance
6. The participants will have access to the study team if they have any issues with the phone or application
7. If the participant does not answer the phone or come to the 2-week visit to return the phone:
 - a. Call the participant with any alternative phone numbers (listed on the locator form completed at the enrollment visit)
 - b. Log in with the study site google account to: <https://www.google.com/android/find>
 - i. This will provide info on the phone's location and if needed allows you to deactivate the phone
 - c. Perform a home visit as feasible
 - d. If the participant is unable to be contacted and is lost to follow up
 - i. If the phone is loaded with data, it should continuously be uploading while on
 - ii. If the phone or participant is lost and is unlikely to be found, remotely deactivate the phone as above

N. Frequently Asked Questions

1. Can you have multiple phones on the same google account?

Yes, you can have multiple devices with the same account

2. How near should the phone be to the participant?

The phone should ideally be on the participant at all the times (within 1-1.5 meters), such as with a carrying bag or an arm band. At night, they can keep the phone next to them. **Can participant take the phone with them to public or my home toilets?**

Yes. The phone app only uses cough or lung sounds and not any other sounds.

4. When can the phone be more than 1.5 meters away from participant?

When participant takes a bath or shower, participant should take off phone and store in a secure, dry place which can be more than 1.5 meters away. When bathing is finished, participant should put phone on participant's person or keep within 1.5 meters.

5. How long should the participant keep the phone near them for recordings?

Participants should keep the phone on, charged and have the app recording for two weeks. Participants should return the phone in two weeks. Study staff will contact participants to retrieve phone

6. What to do if participant phone was lost, damaged, stolen or no longer works properly?

Participants are encouraged to contact the study team to report any problems with the phone without fear of being blamed for event.

7. Can the participant take the phone to work or wear it in living area around other people in the house?

Yes. Participants are encouraged to wear the device during their normal daily activities. The audio recordings will be specific for cough sound recordings only and not record any conversations.

8. Should I still carry the phone when I wear a face mask for COVID or TB protection?

Yes. Participants should keep the phone with active recording on their person when in public while wearing face masks protective against TB or COVID. Participants should notify study staff if wearing a mask for the majority of daily work so staff can make a note of this.

9. There is someone else at participant work or living area that also has a cough. Will the recordings be mixed?

Participants should be mindful of coughs from other people. Participants should be encouraged to keep a zone of 1.5 meters away from anyone else who coughs regularly at home or in workplace. If participant partner coughs at night, phone should be positioned within 1.5 meters of participant, but phone positioned further than 1.5 away from partner.