

S1: DETAILED EXTRACTION PROTOCOL

Reagents

- Ampure XP paramagnetic beads 450 ml, Beckman-Coulter, #A63882
- Buffer RLT 220 ml, Qiagen, #79216
- Buffer EB 250 ml, Qiagen, #19086, or 10 mM Tris-EDTA pH 8 RNase-free, Ambion/ThermoFisher #AM9849
- Ethanol, absolute, molecular biology grade, Fisher Scientific, #BP2818-4
- Nuclease-free water not DEPC treated, Ambion/ThermoFisher, #AM9937
- Bleach, 100%

Consumables

- Abgene™ 96 well 0.8 mL Polypropylene Deep-well Storage Plate - # AB0765
- Eppendorf twin.tec® PCR Plate 96, skirted, 250 µl, PCR clean, colorless, 25 plates – 0030.129.512
- Eppendorf Combitips advanced, 10 mL, 0030.089.820
- Axygen Reagent reservoir, RES-V-25-SI
- VWR Falcon tube 50 mL

Instruments

- Biosafety level 2 cabinet
- Multichannel pipettes range 1-10 µl, 10-200 µl, 100-1000 µl, or similar
- Magnetic plate: Alpaqua #A001322 96S Super Magnet Plate, spring-loaded, or DynaMag-96 side skirted magnetic plate, Life Technologies, #12027
- QuantStudio 12K Flex Real-Time PCR System (OR) Applied Biosystems® 7500 Fast Real-Time PCR System
- HAMILTON Microlab STAR Liquid Handling System with 8-Channel CO-RE (Compressed O-Ring Expansion) Head, 96 well plate Heater Shaker & DynaMag-96 side skirted magnetic plate, Life Technologies, #12027.

Protocol

1. Prepare a fresh 80% Ethanol solution from Ethanol absolute and nuclease-free water.
2. Make sure AmpureXP beads are equilibrated to room temperature. Vortex or mix the beads well before use!
3. For each sample to be extracted, add 150 µl of buffer RLT with 0.5 µg Carrier RNA to the assay plate (Abgene™ 96 Well 0.8 mL Polypropylene Deep-well Storage Plate)
4. Transfer the assay plate to BSL2 cabinet for adding COVID-19 samples.
5. Add 150 µl of sample VTM (Viral Transport Medium) to the assay plate and mix well by pipetting the entire volume up and down 5 times and leave the plate at RT for 5 min to allow lysis.
6. Remove the assay plate from the BSL2 cabinet. Change your gloves and PPE, wipe the outside of the plate with disinfectant to make sure surfaces are clean.
7. Further processing of assay plate can be done on a Hamilton STAR Liquid Handling System (auto) or manually (man).
 - a. (auto) The assay plate is transported to heater shaker for circular shaking at 1600 rpm for 2 min. After shaking, the assay plate is transported back to deck and incubated for 4 min at RT.
 - b. (man) mix the sample by manually pipetting up and down for 10 times and incubate at RT for 4 min.
8. To each sample, add 300 µl of AmpureXP beads and
 - a. (man) mix well by slowly pipetting up and down for 10 times and incubate the assay plate for 9 min or
 - b. (auto) mix cycles - 10, volume - 350 µl, height 2 mm.

9. Transfer the assay plate onto magnetic plate and magnetize the samples for 8-10 min or until the solution appears clear. The exact time depends on the specific plate type and magnet used.
10. Remove and discard the supernatant (595 μ l) without disturbing the bead pellet.
11. Add 300 μ l of freshly prepared 80% Ethanol to each well and mix the pellet by
 - a. (auto) shaking it at 1800 rpm for 120 s on Hamilton Heater Shaker.
 - b. (man) slowly pipetting it up and down until the pellet is fully resuspended.
12. Transport the assay plate back to the magnet and incubate the samples for 6 min or until the solution appears clear.
13. Remove and discard the supernatant without disturbing the bead pellet.
14. Add 300 μ l of freshly prepared 80% Ethanol to each well. DO NOT MIX! Incubate for 30 s. Remove and discard the supernatant without disturbing the bead pellet. Verify the presence of residual ethanol!
15. Remove the plate from the magnet and air dry the pellet for 7 min.
16. Add 23ul of RSB to each well and
 - a. (Auto)resuspend the pellet by shaking it at 1800 rpm for 120 s on Hamilton Heater Shaker or
 - b. (man) pipetting the sample up and down until the pellet is dissolved.
17. Transport the assay plate back to deck and incubate the plate for 4min.
18. After incubation Transport the assay plate back to the magnet and incubate the samples for 3-5 min or until the solution appears clear.
19. Aspirate 20 μ l of the supernatant to a new PCR plate for downstream applications.