


Self-collect nasal and throat swab instructions for COVID-19


Nasal swab

- 1) Wash hands.
- 2) Ensure the swab is labelled with surname, first name, date of birth, gender, time and date of collection, and **site of collection**.
- 3) Gently unscrew the orange lid and remove the swab from the container.
- 4) Do not touch any part of the swab other than the top of the swab.
- 5) Insert the tip of the swab 2-3 cm horizontally into one nostril, as is comfortably possible.
- 6) Rotate the swab **5 times and then leave in position for 5-10 seconds** to absorb any material present.
- 7) Withdraw the swab carefully, taking care not to touch the exposed tip.
- 8) Insert swab shaft into tube, snap shaft at marked line and recap tube.
- 9) Screw on lid firmly to ensure no leakage.



Throat swab

- 1) Wash hands.
- 2) Ensure the swab is labelled with surname, first name, date of birth, gender, time and date of collection, and **site of collection**.
- 3) Gently unscrew the orange lid and remove the swab from the container.
- 4) Do not touch any part of the swab other than the top of the swab.
- 5) In front of a mirror, tilt your head back slightly and stick your tongue out.
- 6) If necessary, use a tongue depressor or spoon to push the back of the tongue down and expose the tonsil area.
- 7) Without touching the sides of the mouth, swab the back of the throat and the tonsil area on both sides of the throat.
- 8) Withdraw the swab carefully, taking care not to touch any other part of the mouth.
- 9) Insert swab shaft into tube, snap shaft at marked line and recap tube.
- 10) Screw on lid firmly to ensure no leakage.



Packaging completed swabs for transport

- Place both labelled swabs into the biohazard bag provided and seal the opening using the ziplock feature.
- Ensure all your personal details on the pathology request form are correct, then fold the form and place into the side pocket on the biohazard bag.

Delivery instructions

Please return your swabs for testing by placing them in the box marked 'SELF-COLLECT SWABS DROP-OFF' located in any of the dedicated COVID-19 collection centres listed on our website, www.dhm.com.au





COVID-19 virus swab collection

Kit contents

Both a nasal and throat swab are required for optimum testing for COVID-19 virus. Instructions on how to collect these specimens are provided overleaf.

- 1 x Nasal swab
- 1 x Throat swab
- 1 x Tongue depressor
- 1 x Specimen bag
- 1 x Instructions

Please visit your
GP to obtain a valid
pathology request
form for COVID-19
testing.

