

Supplementary Material

Material and Methods

RNA Extraction

RNA extraction was performed using a custom automated method for the Analytik Jena CyBio FeliX using the Maxwell HT Viral TNA Kit, Custom (Promega). Sample volumes used were either 500 µl of inactivated sample directly from a Cobas® PCR tube or 200 µl of sample in viral transport media inactivated with 300 µl of the Maxwell Viral Lysis Buffer.

Library Preparation for Illumina Sequencing

Library preparation was performed using a protocol adapted from the COVID-19 ARTIC v3 Illumina Library Construction and Sequencing Protocol v2 (available <https://www.protocols.io/groups/coguk/publications>). In brief, cDNA was synthesised from extracted RNA of combined oropharyngeal and nasopharyngeal swab samples from 16 COVID-19 positive residents and one staff member (tested by RT-qPCR), using LunaScript RT SuperMix (NEB). Enrichment of SARS-CoV-2 DNA was performed using the ARTIC Network nCoV-2019 v3 primers (IDT). Amplified cDNA was purified using 0.8X volume AMPure XP beads (Beckman Coulter) and quantified using the Qubit dsDNA HS Kit (ThermoFisher Scientific). Libraries were prepared using the NEBNext Ultra II DNA Library Prep Kit for Illumina (NEB) and xGen Stubby Adapter and UDI Primer Pairs (IDT). Individual libraries were purified with 0.8X volume AMPure XP beads (Beckman Coulter) and quantified using the Qubit dsDNA HS Kit before constructing normalised library pools. Normalised pools were purified and quantified again before sequencing on a MiSeq (v2 kit, 500 cycles) with 10% PhiX.

Sequence Analysis

Illumina MiSeq paired-end sequencing reads were assembled against the Wuhan-Hu-1-2019 reference genome (NCBI accession MN908947) using BWA-MEM (Li 2013, <https://arxiv.org/abs/1303.3997>). Single nucleotide variants (SNV) were called in the remaining 17 sample sequences using samtools with minimum depth of 100 reads. Unrooted phylogenetic trees were generated from SNV data using PHYLIP (dnadist F84 and neighbour-joining).